



7 Leonard Street, Hornsby NSW 2077 Australia-P.O. Box 1664, Hornsby Westfield NSW 1635 Tel: (02) 9476 4544 • Fax: (02) 9477 7974 • Int'l +61-2 • E-mail: kelba@kelba.com • Web Site: www.kelba.com

Kelba KH-MH500 Wheel Chair Platform – Scale

Kelba's KH-MH500 wheel chair scale offers the user easy and safe access via built-in entry and exit ramps, reducing patient disturbance. Portable design includes wheels and handles located on its side to make moving the scale from place to place an easy task. Ideally suited in a vast range of medical applications for either use with wheel chairs or as a stand on scale. Its 500kg capacity can cover a wide range of weighing applications. The scale includes either a pole or desk / wall mountable indicator with an 25mm LCD display, Backlight, with low power and over load indication.







7 Leonard Street, Hornsby NSW 2077 Australia-P.O. Box 1664, Hornsby Westfield NSW 1635 Tel: (02) 9476 4544 • Fax: (02) 9477 7974 • Int'l +61-2 • E-mail: kelba@kelba.com • Web Site: www.kelba.com

Kelba KH-MH500 Wheel Chair Scale

Key Features

For wheel chair and stand on weighing	Built-in ramps
Extra wide for large wheel chairs	Pole mounted indicator
Wheels & handles for portability of scale	🔸 🛛 Backlit LCD
🖕 Multiple uses	Hold Function
Aluminium Build	RS232 interface

Specifications

Model	Capacity	Increments	Display	Power Requirements
КН-МН500	500kg	200 grams	LCD 25mm	AC power &
	300kg	100 grams	Backlit	DC 6 volt rechargeable

Finish:	Non-Slip Aluminium
Operating Temperature:	5°C to 40°C
Storage Temperature:	0°C to 50°C
Operating Humidity:	85%R.H
Display Reading:	Kg, g, lb
Physical Foot Print Dimensions:	1200mm long x 1050mm wide x 100mm High with built-in ramps
Platform Size:	800mm long x 800mm wide x 50mm High
Net Weight:	31kg
Shipping Weight:	36kg